

# A clinical trial of cResponse, a functional assay for cancer precision medicine

Seth Salpeter<sup>1</sup>, Chani Stossel<sup>4</sup>, Vered Bar<sup>1</sup>, Sara Aharon<sup>1</sup>, Luba Torovsky<sup>1</sup>, Adi Zundelevich<sup>1</sup>, Hamutal Shahar<sup>1</sup>, Hagit Shapira<sup>1</sup>, Giuseppe Mallel<sup>1</sup>, Nancy Gavert<sup>2</sup>, Ravid Straussman<sup>2</sup>, Shay Golan<sup>3</sup>, Eli Rosenbaum<sup>3</sup>, Zohar Dotan<sup>4</sup>, Raanan Berger<sup>4</sup>, Dan Leibovici<sup>5</sup>, Shani Breuer<sup>6</sup>, Yakir Rottenberg<sup>6</sup>, Ayala Hubert<sup>6</sup>, Hovav Nechushtan<sup>6</sup>, Aviad Zick<sup>6</sup>, Guy Neev<sup>1</sup>, Talia Golan<sup>4</sup>


<sup>1</sup>Curesponse, Weizmann 6, Tel Aviv Israel. | <sup>2</sup>Department of Molecular Cell Biology, Weizmann Institute of Science, Rehovot Israel. | <sup>3</sup>Department of Urology, Beilinson Hospital – The Rabin Medical Center, Petach Tikva Israel.

<sup>4</sup>Department of Oncology, Sheba Hospital, Ramat Gan Israel. | <sup>5</sup>Department of Urology, Kaplan Medical Center, Rehovot Israel. | <sup>6</sup>Sharett Institute of Oncology, Hadassah Medical Center, Hebrew University Jerusalem Israel.


## Abstract

Precision cancer therapy has the potential to revolutionize treatment and improve outcome. While genomic analysis has become central to cancer personalized medicine, recent studies have not shown that it drastically improves the patient's survival as compared to drug selection. Additionally, the presence of genomic mutations may suggest several treatment options without elucidating which of the drugs or combination of drugs might yield superior clinical response. Moreover, for many drugs no genetic predictive biomarkers are available. To advance cancer precision diagnostics, we have developed cResponse, a functional drug sensitivity platform to determine individualized patient treatment regimens. Fresh patient cancer tissue samples are obtained by biopsy or resection and sectioned to 300 μm slices which when cultured in the cResponse platform, retain their architectural character and proliferative capacity. We show that cResponse preserves human cancer tissue together with its microenvironment, including endothelial and immune cells, at a high viability (>90%) with continued cell division for 7 days in over twenty types and sub-types of solids tumors. We utilize this extended period of viability to (a) perform rapid genomic profiling to help prioritize drugs to be tested by cResponse and (b) evaluate the effects of slow-acting drugs such as targeted therapy. To test the accuracy of the cResponse assay to predict patient response to cancer treatment, we first conducted a comparative study on thirty-seven patients with metastatic pancreatic cancer receiving platinum-based drugs and Olaparib. The tissues were treated for 5 days ex vivo and assessed with histological staining. Tissues were assigned a viability score (ranging from 0 to 100) based on an algorithm composed of a panel of immunohistochemical markers. Clinical correlation showed a specificity of 90% and a sensitivity of 88%. Subsequently, we performed a clinical study on biopsies from twenty-eight patients with different cancer types receiving diverse chemotherapeutic and targeted agents. Cancer types included bladder, pancreatic, lung, colorectal, breast and sarcoma. The cResponse score was compared to the clinical response and showed that cResponse could predict the patient's response with a specificity of 78% (7/9) and a sensitivity of 95% (18/19). **Taken together, a total of 65 samples demonstrated a specificity of 87% (33/38) and a sensitivity of 96% (25/27).** In the future, the integration of this platform in the clinical setting, to help direct drug choices for anti-cancer treatment may lead to improved patient outcomes.


## cResponse Ex Vivo Organ Culture System Preserves the Tumor Microenvironment



(A) Resected Colorectal Cancer (CRC), Transitional Cell Carcinoma (TCC) and Breast Cancer (BC) were cultured in an ex vivo organ culture assay capable of preserving the tumor microenvironment. Representative images of the tissues are compared on Day 0 when the tissue was received and after 5 days of culture. H&E and Ki67 stains are shown.




(B) Viability and proliferation of the different cancer samples were compared between Day 0 and Day 5 (N=4). Cancer viability and ki67% stain were assessed and quantified by a pathologist showing non-significant changes from the beginning of the assay.



(C) The effectivity of targeted drug treatments to modulate their targets in ex vivo culture was assessed after 24 hours. Tissues were treated with either Trametinib, a pERK inhibitor, Palbociclib, a CDK4/6 inhibitor or NT219, a pStat3/IRS inhibitor as indicated. Their respective molecular targets were evaluated by IHC staining.

## Comparative Study




(A) A comparative study was conducted using cResponse to predict the clinical response of patients with metastatic pancreatic cancer. Tumor section from PDAC PDX models were assessed by cResponse and correlated to the original patient's clinical response to platinum and PARPi at the time of PDX generation. Representative images are shown of a tissue sample after 5 days in cResponse with vehicle and cisplatin.

Figure B is a summary table of cResponse scores compared to clinical response for 37 samples. A central callout states: "cResponse score was compared to the clinical response for 37 samples showing: 90% specificity (26/29) and 88% sensitivity (7/8)".


Patient	Score	Clinical Response	Treatment
126b	0	No-Response	Olaparib
217b	0	No-Response	Olaparib
321a	4.5	No-Response	Cisplatin
187b	7	No-Response	Olaparib
126a	8.5	No-Response	Cisplatin
187a	9	No-Response	Cisplatin
418a	14	No-Response	Cisplatin
385b	16	No-Response	Olaparib
291a	21.5	No-Response	Cisplatin
291b	22.9	No-Response	Olaparib
227b	24	No-Response	Olaparib
194a	24	No-Response	Cisplatin
210a	24.7	No-Response	Cisplatin
143b	26.7	No-Response	Olaparib
375a	29	No-Response	Cisplatin
185a	32	No-Response	Cisplatin
421a	32	No-Response	Cisplatin
210b	33.7	No-Response	Olaparib
421b	37	No-Response	Olaparib
265a	49	No-Response	Cisplatin
265b	49.1	No-Response	Olaparib
143a	53.7	No-Response	Cisplatin
217a	61	No-Response	Olaparib
145b	61.1	No-Response	Olaparib
96a	63.8	Response	Olaparib
135a	66	No-Response	Cisplatin
228a	66	No-Response	Cisplatin
145a	66.5	Response	Cisplatin
326a	69	Response	Olaparib
259a	69	Response	Cisplatin
228b	74	No-Response	Olaparib
203a	76	Response	Cisplatin
203b	82	Response	Olaparib
122a	87	Response	Cisplatin
87a	87	No-Response	Olaparib
324a	90	No-Response	Cisplatin
467a	100	Response	Cisplatin

(B) Summary of cResponse score compared to the clinical response based on oncologist evaluation (sensitive - SD, PR, CR; resistance - PD). Patients received platinum-based therapy followed by PARPi. Response to each treatment was recorded independently. cResponse score was assessed on PDX tumor chunks recapitulating the specific clinical scenario and response to platinum and PARPi.

## Clinical Study



(A) Biopsies from breast cancer (BC), pancreatic cancer (PDAC) and sarcoma were obtained prior to initiation of patient treatment and maintained in the ex vivo organ culture for 5 days. Representative images are shown of each cancer type at Day 0 and Day 5 showing high viability and preservation during the culture assay.



(B) Representative images of samples from the clinical trial which were responders and non-responders. The responder is from a pancreatic cancer patient treated with Gemzar-Taxol, while the non-responder is from a bladder cancer patient treated with Cisplatin and Gemzar.

Figure C is a table of results from the cResponse test and clinical correlation. A central callout states: "cResponse score was compared to the clinical response for 28 samples showing: 78% specificity (7/9) and 95% sensitivity (18/19)".

Patient	Score	Clinical Response	Tissue Type	Treatment
e132	0	PD	Bladder	Cis-Gemzar
55	8	PD	Bladder	Cis-Gemzar
25	8	PD	Bladder	Cis-Gemzar
e55	13	PD	Sarcoma	Paoparib-Everolimus
c229	13	PD	Esophageal	Paclitaxol
70	17	PD	Bladder	Cis-Gemzar
c150	31	PR	CRC	Folfox-Cetuximab
e93	37	PD	Pancreatic	FOLFIRINOX
C139	45	SD	Bladder	Cis-Gemzar
159	46	SD	Sarcoma	Ifosamide-Etoposide
c161	50	SD	Pancreatic	FOLFIRINOX
117	62	SD	Pancreatic	FOLFIRINOX
156	64	PR	Pancreatic	Gemzar-Taxol
107	66	SD	Pancreatic	FOLFIRINOX
130	74	PD	Bladder	Cis-Gemzar
e83	75	PR	Bladder	Cis-Gemzar
170	82	PD	Unknown	Cis-Gemzar
41	86	PR	Bladder	Cis-Gemzar
149	86	PR	CRC	FOLFOX
C168	90	PR	Pancreatic	FOLFIRINOX
115	91	PR	Pancreatic	FOLFIRINOX
128	100	PR	Bladder	Cis-Gemzar
190	92	PR	Bladder	Cis-Gemzar
33	100	PR	Bladder	Cis-Gemzar
63	100	PR	Bladder	Cis-Gemzar
95	100	PR	Bladder	Cis-Gemzar
105	100	PR	Bladder	Cis-Gemzar
c135	100	CR	Breast	ACT-C

(C) Table of the results from the cResponse test and clinical correlation, leading to overall sensitivity and specificity calculations. Thresholds for response were 0-40 for progressive disease, 40-80 for stable disease & partial response and 80-100 for partial response and complete response. Clinical response was determined by either pathological or clinical imaging.

## Conclusions

- cResponse assay maintains high viability of the resected human cancer tissue, preserving the 3D architecture together with its microenvironment, including endothelial and immune cells for 5 days.
- Cell signaling pathways within the cancer tissue are maintained as demonstrated by the downregulation of target proteins after treatment ex vivo.
- Similar to resected tissue, we show the ability to preserve core biopsies at high viability in our system and evaluate their response to drugs.
- In a comparative study in pancreatic metastatic cancer, the cResponse score was compared to the clinical response for 37 samples, showing 90% specificity and 88% sensitivity.
- An expanded clinical study was performed on biopsies from patients with locally advanced and metastatic cancer in patients receiving diverse chemotherapeutic and targeted anti-cancer agents. Cancer types included bladder, pancreatic, lung, colorectal, breast and sarcoma. In 28 biopsy samples, the cResponse assay demonstrated 78% specificity and 95% sensitivity when cResponse assay as compared to clinical outcome.

A total of 65 samples demonstrates specificity of 87% (33/38) and sensitivity of 96% (25/27)

## References

- Ben-Hamo R, Jacob Berger A, Gavert N, et al. Predicting and affecting response to cancer therapy based on pathway-level biomarkers. Nat Commun. 07 2020;11(1):3296.doi:10.1038/s41467-020-17090-y
- Nejman D, Livyatan I, Fuks G, et al. The human tumor microbiome is composed of tumor type-specific intracellular bacteria. Science. 05 2020;368(6494):973-980. doi:10.1126/science.aay9189
- Grinshpun A, Gavert N, Granit RZ, et al. Ex vivo organ culture as potential prioritization tool for breast cancer targeted therapy. Cancer Biol Ther. Mar 2018;1-4.doi:10.1080/15384047.2018.1450114